

# Data Sheet



## Excel Life™ Medical Range Common Plate, Double Powerpoint, Antimicrobial\*, 10A

marked "RCD PROTECTED"

Cat No. **EC777CAWE, EC777CABU, EC777CARD**



<b>Cat No.</b>	<b>EC777CAWE, EC777CABU, EC777CARD</b>
<b>Description</b>	Excel Life Medical Range - Common Plate, Double Powerpoint, Antimicrobial, 10A
<b>Supply voltage</b>	230 – 240V a.c.
<b>Frequency</b>	50Hz
<b>Maximum current</b>	10A
<b>Colours</b>	White (WE), Blue (BU), Red (RD)
<b>Dimensions (WxHxD)</b>	118 x 74 x 10mm
<b>Mounting centres</b>	84mm
<b>Operating temperature</b>	5° to 35°C
<b>Cut-out dimensions</b>	56 x 96mm
<b>Standard</b>	AS/NZS 3112

# Excel Life™ Medical Range Common Plate, Double Powerpoint, Antimicrobial\*, 10A marked "RCD PROTECTED"



Cat No. **EC777CAWE, EC777CABU, EC777CARD**



- Chemical Resistant and Antimicrobial
- Internal safety shutters for extra safety
- More fixing options – 4 corners
- Captive terminal screws
- Marked "RCD Protected"

*Antimicrobial Specifications			
<b>Maintenance</b>	Surface cleaning with cloth: Hospital-grade water-based cleaners only. Do not use any solvent based cleaner such as acetone, tar-remover or trichloroethylene.		
<b>Antimicrobial Properties</b>	Challenge Organism	Polycarbonate alloy in facia parts	Polycarbonate alloy in switch dolly
		Reduction Coefficient/Kill Rate	
<b>Note:</b> 5 Log = 99.9999% 4 Log = 99.99% 3 Log = 99.9% 2 Log = 99%	<b>Escherichia Coli</b>	2.55 Log / 99.72%	6.72 Log / 99.9999%
	<b>Staphylococcus aureus</b>	2.70 Log / 99.79%	5.22 Log / 99.9999%
<p>Tests performed according to JIS Z 2801:2010 (Antimicrobial products – Test for antimicrobial activity and efficacy). Test method is equivalent to procedure detailed in ISO 22196:2011 (Measurement of antibacterial activity on plastics and other non-porous surfaces). Test results available upon request.</p> <p>Not to be used as an alternative or substitute for standard hygienic practices, procedures, or protocols.</p> <p>This device has a long-term antimicrobial effect. Lifespan is dependent upon cleaning procedures and the use of recommended cleaning products.</p> <p>Product manufactured with antimicrobial material based on Ag+ ions. This technology limits the proliferation of bacteria on our devices without increasing antibiotic resistance.</p>			

Single product	
Cat No.	Barcode
EC777CAWE	9321001428251
EC777CABU	9321001442448
EC777CARD	9321001442455

Inner box (Qty 10)		
Cat No.	Barcode	Dimensions (WxHxD)
EC777CAWE	9321001429067	175 x 165 x 125mm
EC777CABU	9321001443193	175 x 165 x 125mm
EC777CARD	9321001443209	175 x 165 x 125mm

Outer box (Qty 120)		
Cat No.	Barcode	Dimensions (WxHxD)
EC777CAWE	9321001429876	608 x 308 x 340mm
EC777CABU	9321001443919	608 x 308 x 340mm
EC777CARD	9321001443926	608 x 308 x 340mm

Legrand - Australia  
1300 369 777  
www.legrand.com.au  
ABN 31 000 102 661

Legrand - New Zealand  
0800 476 009  
www.legrand.co.nz

Refer to your group buying office, Legrand Sales Representative or Legrand Sales Office for pricing.

\* Specifications may change without notice