

Data Sheet



Excel Life™ Medical Range

Common Plate, Double Powerpoint
with Extra Function Hole, Antimicrobial*, 10A
marked "RCD PROTECTED"

Cat No. **EC777XCAWE, EC777XCABU, EC777XCARD**



Cat No.	EC777XCAWE, EC777XCABU, EC777XCARD
Description	Excel Life Medical Range - Common Plate, Double Powerpoint with Extra Function Hole, Antimicrobial, 10A
Supply voltage	230 – 240V a.c.
Frequency	50Hz
Maximum Current	10A
Colours	White (WE), Blue (BU), Red (RD)
Dimensions (WxHxD)	118 x 74 x 10mm
Mounting centres	84mm
Operating temperature	5° to 35°C
Cut-out dimensions	56 x 96mm
Standard	AS/NZS 3112

Excel Life™ Medical Range Common Plate, Double Powerpoint with Extra Function Hole, Antimicrobial*, 10A marked "RCD PROTECTED"



Cat No. **EC777XCAWE, EC777XCABU, EC777XCARD**



- Chemical Resistant and Antimicrobial
- Internal safety shutters for extra safety
- More fixing options – 4 corners
- Extra function hole

*Antimicrobial Specifications			
Maintenance	Surface cleaning with cloth: Hospital-grade water-based cleaners only. Do not use any solvent based cleaner such as acetone, tar-remover or trichloroethylene.		
Antimicrobial Properties	Challenge Organism	Polycarbonate alloy in facia parts	Polycarbonate alloy in switch dolly
		Reduction Coefficient/Kill Rate	
Note: 5 Log = 99.999% 4 Log = 99.99% 3 Log = 99.9% 2 Log = 99%	Escherichia Coli	2.55 Log / 99.72%	6.72 Log / 99.9999%
	Staphylococcus aureus	2.70 Log / 99.79%	5.22 Log / 99.9999%
	Tests performed according to JIS Z 2801:2010 (Antimicrobial products – Test for antimicrobial activity and efficacy). Test method is equivalent to procedure detailed in ISO 22196:2011 (Measurement of antibacterial activity on plastics and other non-porous surfaces). Test results available upon request. Not to be used as an alternative or substitute for standard hygienic practices, procedures, or protocols. This device has a long-term antimicrobial effect. Lifespan is dependent upon cleaning procedures and the use of recommended cleaning products. Product manufactured with antimicrobial material based on Ag+ ions. This technology limits the proliferation of bacteria on our devices without increasing antibiotic resistance.		

Single product	
Cat No.	Barcode
EC777XCAWE	9321001428268
EC777XCABU	9321001442448
EC777XCARD	9321001442455

Inner box (Qty 10)		
Cat No.	Barcode	Dimensions (WxHxD)
EC777XCAWE	9321001428268	85 x 117 x 50mm
EC777XCABU	9321001442448	85 x 117 x 50mm
EC777XCARD	9321001442455	85 x 117 x 50mm

Outer box (Qty 120)		
Cat No.	Barcode	Dimensions (WxHxD)
EC777XCAWE	9321001429883	608 x 308 x 340mm
EC777XCABU	9321001443919	608 x 308 x 340mm
EC777XCARD	9321001443926	608 x 308 x 340mm

Legrand - Australia
1300 369 777
www.legrand.com.au
ABN 31 000 102 661

Legrand - New Zealand
0800 476 009
www.legrand.co.nz

Refer to your group buying office, Legrand Sales Representative or Legrand Sales Office for pricing.

* Specifications may change without notice